Publication

UV VIS

Fields of Application / Industry:

- Chemistry / Polymer Industry
- Clinical Chemistry / Medicine / Hygiene / Health Care
- Cosmetics
- Electronics
- Energy
- Environment / Water / Waste
- Food / Agriculture
- Geology / Mining
- Material Analysis
- Metallurgy / Galvanization
- Pharmacy
- Refineries / Petrochemistry
- Semi-Conductor Technology
- Others

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Enzymes – Biological Catalysts in Food Analysis

Alexandra Kästner, Application Chemist Molecule Spectroscopy, Analytik Jena AG,
Analytical Instrumentation, Konrad-Zuse-Str. 1, 07745 Jena

Introduction

As biological catalysts, enzymes control all metabolic processes in the human organism. They are, among other things, responsible for the digestion of nutrients, the regeneration of cells and tissue, the function of individual organ systems and for the biosynthesis of proteins as well as hormones. They support the immune system and are involved in detoxification processes. Without enzymes all metabolic processes would be impossible. From a chemical point of view they are protein molecules with functional groups which interact with the substance to be analyzed (the substrate) and on which the substance conversion is performed. During this, the substrate is either converted itself via the enzyme or it influences the enzyme-catalyzed reaction. Enzymes have a distinctive selectivity for binding certain molecules and always catalyze a defined reaction, e.g. hydrolysis, oxidation/reduction or rearrangement, which can be deviated from its name (e.g. hydrolase, oxidase). [1] Because of these characteristics, enzymatic methods offer significant benefits compared to the chemical processes when examining countless materials such as food, biological samples and many others. They are used in different areas of the industry, research and official monitoring.

In food chemistry, enzymatic processes are used to evaluate foods. Ingredients, such as sugar, acids and their salts, alcohol or other substances, can quickly be determined quantitatively even from complex matrices by performing the relevant sample preparation. Furthermore measurements of the activity of certain enzymes give an indication of microbial pollution or can be used to test a sufficient heating of foods during pasteurization or sterilization processes because of their sensitivity to heat [2].



Experimental

In the following application example the concentration of sulfite in white wine is determined using enzymatic determination. Sulfite or sulfur dioxide is e.g. added to wine as antioxidant and preservative or to bind fermentation by-products which would have a negative impact on the taste. The total amount of sulfur dioxide (total of bound and free SO_2) that wine may contain is limited by

an EU regulation and depends on the type, the quality level as well as the sugar content of the wine. [3]

For the measurement the SPECORD® 210 PLUS with a Peltier-tempered 8-cell changer with external heat exchanger was used. (Fig. 1) The SPECORD® 210 PLUS is a double-beam photometer for the wavelength range of 190–1100 nm, with 4-fold variable spectral slit width and two tempered photo diodes (CDD – Cooled Double Detection). [4] The Peltier-tempered 8-cell changer is an autosampler with 8 positions for cells with a path length of 10 mm (Fig. 2). It enables enzymatic analyses of multiple samples which require very exact tempering. The cell block is connected to an external heat exchanger via hoses, which is used to cool the back of the Peltier element. The temperature control is performed by a separate control unit, for which the temperature can be controlled in a range between -5°C and +105°C with an accuracy of +0.1°C. A measurement sensor is used as control sensor and positioned at the bottom outer corner of the cell block. The cell changer further features two additional sensors for the registration of the holder or inside cell temperature respectively. A measurement sensor is used for the direct control of the inside cell temperature which can remain in the cell during the optical measurement. [5]

The enzymatic reactions are performed in several steps. In the first step of the reaction the sulfite (SO_3^{2-}) is oxidized in the presence of oxygen by the enzyme sulfite oxidase ($\text{SO}_2\text{-OD}$) to sulfate (SO_4^{2-}) (1). Hydrogen peroxide (H_2O_2) which is created in this step is converted in the second step with reduced NADPH by the enzyme NADPH peroxidase (NADPH-POD) to water and NAD^+ (2). The amount of used NADH is measured photometrically at 340 nm and is equivalent to the sulfite concentration. [6]

A test combination of R-Biopharm [6] and the sulfite method of the WinASPECT PLUS® software were used for the determination of the concentration of sulfite, which is performed at a temperature of 25°C. The white wine did not require any further sample preparation and could be used with 0.1 ml directly with the reaction solutions for the measurement. According to the instruction of the test kit 1.0 ml of the reaction mixture 2 (consisting of triethanolamine buffer and NADPH) was placed into the cells, mixed for the blank value with 2.0ml distilled water, for the sample with 1.9 ml distilled water and then 0.01 ml of suspension 3 (consisting of NADH-POD) were added. The reference measurement was performed with water. All reagents were mixed with a stirring spatula

and after an incubation period of 5min the absorbance E_1 could be measured. The enzymatic reaction was started by adding 0.05 ml of the suspension 4 (SO_2 -OD). After a further incubation period of 30 min the second absorbance (E_2) was measured. Additionally a sulfite control standard was measured to check the accuracy of the analysis. The processing of the total method is performed automatically by the WinASPECT PLUS[®] software.

Results and evaluation

The absorbance differences for the blank value and the samples are calculated to determine the concentration of sulfite. Afterwards the absorbance difference of the blank value is subtracted from the absorbance differences of the samples:

$$\Delta E = (E_1 - E_2)_{\text{Sample}} - (E_1 - E_2)_{\text{Blank value}}$$

The sulfite concentration is calculated according to the general formula:

$$c \text{ [g/l]} = \frac{V \cdot MG}{\varepsilon \cdot d \cdot v \cdot 1000} \cdot \Delta E$$

whereby	V	= Test volume [ml]
	v	= Sample volume [ml]
	MG	= Molecular weight of the substance to be determined [g/mol] SO_2 64.06 g/mol
	d	= Path length of the cell [cm]
	ε	= Absorption coefficient of NADH at 340 nm = $6.3 \text{ [l} \cdot \text{mmol}^{-1} \cdot \text{cm}^{-1}]$

The calculation of the absorbance differences and the sulfite concentration are performed automatically by the software after the measurement. A concentration of sulfite of 60.1 mg/l was determined for the white wine sample. The maximum quantity of sulfur dioxide in white wine with less than 5 g residual sugar/l is 200 mg/l. According to the Food Labeling Directive of the EU wines which contain a concentration of more than 10 mg/l SO_2 have to be labeled because sulfite belongs to the allergenic substances. Any wine naturally includes SO_2 which is produced in low concentrations (up to 40 mg/l) under normal, physiological conditions during fermentation. [3]

Summary

Because of their effect and substance specificity the enzymes enable a number of simultaneous metabolic processes in the human organism. In the routine analysis enzymatic methods are therefore of great importance. By the measurement of the absorbance change in the UV range of absorbent coenzymes nearly all substances which can be converted by enzymes can be determined photometrically.

An enzyme-based analysis is fast and can be performed without time-consuming sample preparation. It is performed under physiological conditions. The reagents are harmless and easy to handle. Furthermore the documentation of individual substances in mixtures is possible. Thanks to an extensive collection of enzymatic methods in the WinASPECT PLUS® software, the determination of the corresponding substance with the relevant measurement setting, sample sequence and evaluation is performed automatically.

SPECORD® PLUS accessories:

Fig. 1 SPECORD® PLUS with Peltier-tempered cell holder with external heat exchanger


Fig. 2 Peltier-tempered 8-cell changer

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Printout and further use permitted with reference to the source.

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© 2012 Analytik Jena AG